


Patient-Specific PEEK Implants



- Implantable PEEK
- Strong and Durable
- Customized Fit
- Fast Delivery

Surgeons choose **ACCUSHAPE® PEEK Cranial Implants** because of fast delivery, PEEK material, and customized fit.

Uses

- Cranioplasty
- Cosmesis
- Bone Flap Replacement
- Volume / Void Filler

Benefits


- **Customizable** design
- **Save time** with patient-specific implant designed to fit without modification in surgery
- **Compatible** with titanium neuro plating systems
- **Permanent** alternative to bone flap

PEEK Material

- Strength comparable to cortical bone
- Biocompatible
- Radiolucent

ACCUSHAPE® PEEK Cranial Implants

START	SUBMIT SERVICE REQUEST and CT SCAN(S) , TO MEDCAD
PLAN	PLAN SURGERY VIA ONLINE MEETING (as needed)
	MEDCAD SENDS THE PLAN FOR APPROVAL
APPROVE	SURGEON APPROVES IMPLANT DESIGN or REQUESTS CHANGES
MAKE	MEDCAD MANUFACTURES & SHIPS
GET	GET ACCUSHAPE® PEEK CRANIAL IMPLANT
SURGERY	



Customizable Features

1. Drainage, Perforation, & Tack-Up Holes
2. Pre-Drilled Fixation Pilot Holes
3. Edge Retraction & Contour

For Service Request:
MedCAD.net/start

For Digital File Transfer:
MedCAD.net/upload



Ship To: 501 S 2nd Ave, Suite A-1000, Dallas, TX 75226

Phone: +1 (214) 453-8864

Email: orders@medcad.net

Website: MedCAD.net

